

<b>DATE:</b>	15 September 2025
<b>TO:</b>	Cross Cancer Institute Healthcare Providers
<b>FROM:</b>	Clinical Biochemistry, Alberta Precision Laboratories (APL)
<b>RE:</b>	<b>Implementation of High Sensitivity Troponin I (hs-TnI) and BNP on the Quidel Triage Meter</b>

---

## PLEASE POST OR DISTRIBUTE AS WIDELY AS APPROPRIATE

---

### Key Message

- Effective 16 September 2025, high sensitivity troponin I (hs-TnI) and B-natriuretic peptide (BNP) will be added to the test menu at the Cross Cancer Institute lab on the Quidel Triage Meter. This will improve the turnaround time for these tests.

### hsTnI

- The Triage True hs-TnI assay requires sample collection in lavender EDTA blood collection tubes
  - No other sample collection tubes are acceptable.
- Triage True hs-TnI assay includes a different reference interval (99<sup>th</sup> percentile upper reference limit), reporting limits, rule-in/rule-out chest pain pathway, delta values, critical limits and interpretative comments than the hs-TnT assay used at the University of Alberta Hospital lab (Table 1, Table 2).
- TriageTrue hs-TnI assay will offer the new 2-Hour Chest Pain Pathway for Quidel TriageTrue High Sensitivity Troponin-I Assay (Appendix 1).
- Outside of the Cross Cancer Institute lab operational hours, testing will continue to be performed at the University of Alberta Hospital Lab. Results cannot be compared between methods.
- If hs-TnT is required, this can be ordered with the test orderable: Troponin High Sensitivity Send Out.

### BNP

- The Triage BNP assay requires sample collection in lavender EDTA blood collection tubes.
  - No other sample collection tubes are acceptable.
- Triage BNP assay has different units of measure, reporting limits, reference intervals, and interpretative comments than NT-proBNP assay used at the University of Alberta Hospital lab (Table 3).
- Outside of the Cross Cancer Institute lab operational hours, testing will continue to be performed at the University of Alberta Hospital Lab. Numerical results cannot be compared between methods.

### Background

- Currently these tests are referred out to the University of Alberta Hospital. Addition of these tests to the onsite test menu supports patient care.
- The Ortho Quidel Triage Meter is used across the province in smaller hospital laboratories for troponin and BNP testing.



- Evidence supports that a 2-hour chest pain pathway for TriageTrue hs-Tnl is effective and safe for rule-in/rule-out of acute myocardial infarction (AMI). • The recommended pathway for TriageTrue hs-Tnl is consistent with clinical practice guidelines and is recommended by the Cardiovascular Program Improvement and Integration Network (PIN) in consultation with Emergency Medicine PIN and Laboratory Medicine provincially.

### **How this will impact you**

- Turnaround time will be improved for these tests, supporting patient care at the Cross Cancer Institute.
- There is no change to ordering in ConnectCare.

### **Action Required**

- Ensure sample collection in correct tube type. This will vary depending on Cross Cancer Laboratory operational hours
- Be aware of differences in testing, reporting and interpretation for these tests.

**Effective 16 September 2025**

### **Questions/Concerns**

- Dr. Miranda Brun, Clinical Biochemist, APL: [Miranda.brun2@aplabs.ca](mailto:Miranda.brun2@aplabs.ca); 780-722-4123

### **Approved by**

- Dr. Kareena Schnabl, Section Chief Biochemistry, North Sector
- Dr. Michael Mengel, Medical Director, North Sector



**Table 1: Summary of new reporting changes for troponin at the Cross Cancer Institute**

<b>Testing location</b>	University of Alberta Hospital Lab	Cross Cancer Institute Lab
<b>Test orderable</b>	Troponin	Troponin
<b>Test method</b>	Roche hs-TnT	Quidel TriageTrue hs-TnI
<b>Collection Tube</b>	Barricor lithium heparin PST (lime green)	Lavendar EDTA
<b>Reporting units</b>	ng/L	ng/L
<b>Reference interval (99<sup>th</sup> percentile)</b>	<14 ng/L	<21 ng/L
<b>Critical value</b>	≥ 52 ng/L	≥ 60 ng/L
<b>Reporting limits</b>	3 to 500 000 ng/L	2 to 1000 ng/L
<b>Delta Value</b>	Reported for 0-2 hour delta	Reported for 0-2 hour delta
<b>Rapid Chest Pain Pathway</b>	<a href="#"><u>2-Hour Chest Pain Pathway for Troponin T-High Sensitivity hs-TnT - Roche - Edmonton</u></a>	<a href="#"><u>2-Hour Chest Pain Pathway for Quidel TriageTrue High Sensitivity Troponin-1 Assay</u></a>



**Table 2: Interpretative comments reported with Quidel TriageTrue hs-Tnl**

hs-Tnl result (ng/L)	Comment	Flagging
<4	<p>For patients with a non-ischemic ECG, a Troponin I, High Sensitivity of 3 ng/L or less on presentation is highly sensitive for excluding acute myocardial infarction, provided the specimen was collected more than 3-hours from symptom onset. However, for patients with symptoms less than 3-hours duration or concerning clinical presentations, repeat troponin testing at 2-hours after the initial sample is recommended.</p> <p>Please note that patients with ischemic ECG changes and/or high-risk clinical presentations should be considered for further evaluation irrespective of troponin results.</p>	Normal
4-20	<p>Troponin I, High Sensitivity is below the upper reference limit (21 ng/L) and results are not consistent with myocardial infarction (MI) or injury. However, patients with acute symptoms (less than 6-hours) or concerning clinical presentations should undergo repeat troponin testing at 2-hours after the initial sample.</p> <ul style="list-style-type: none"> <li>- Troponin I, High Sensitivity of 4 ng/L or less on presentation AND a 2-hour delta(change) of 2 ng/L or less is highly sensitive for excluding acute myocardial infarction (MI).</li> <li>- A 2-hour delta (change) of 3-7 ng/L may indicate acute myocardial injury and suggest an additional troponin measurement 4 hours after the initial sample, serial ECG testing and clinical re-evaluation.</li> <li>- A 2-hour delta (change) of 8 ng/L or more suggests an acute myocardial injury and may represent acute myocardial infarction in the appropriate clinical scenario.</li> </ul> <p>Please note that patients with ischemic ECG changes and /or high-risk clinical presentations should be considered for further evaluation irrespective of troponin results</p>	Normal
21-59	<p>Troponin I, High Sensitivity has a non-specific/non-diagnostic elevation. Interpretation is highly dependent on clinical presentation and patient history. New elevations are concerning; however, many patients have chronic elevations in troponin and measured concentrations near the patient's baseline are reassuring. However, patients with acute symptoms (less than 6-hours) or concerning clinical presentations should undergo repeat troponin testing at 2-hours after the initial sample.</p> <ul style="list-style-type: none"> <li>- A 2-hour delta (change) of 2 ng/L or less is highly sensitive for excluding acute myocardial infarction.</li> <li>- A 2-hour delta (change) of 3-7 ng/L may indicate acute myocardial injury and suggest an additional troponin measurement 4 hours after the initial sample, serial ECG testing and clinical re-evaluation.</li> <li>- A 2-hour delta (change) of 8 ng/L or more suggests an acute myocardial injury and may represent acute myocardial infarction in the appropriate clinical scenario.</li> </ul> <p>Please note that patients with ischemic ECG changes and /or high-risk clinical presentations should be considered for further evaluation irrespective of troponin results</p>	High
60	<p>Clear elevation of Troponin I, High Sensitivity consistent with acute myocardial injury or infarction in the appropriate clinical context Repeat troponin testing at 2-hours after the initial sample may be helpful to assess for ongoing myocardial injury</p>	Critical




**Table 3. Summary of new reporting changes for BNP at the Cross Cancer Institute**

<b>Testing location</b>	University of Alberta Hospital Lab	Cross Cancer Institute Lab
<b>Test orderable</b>	B-Natriuretic Peptide (BNP OR NT-proBNP)	B-Natriuretic Peptide (BNP OR NT-proBNP)
<b>Test method</b>	NT-proBNP	BNP
<b>Collection Tube(s)</b>	Mint lithium heparin PST Red Lavendar EDTA Gold SST	Lavendar EDTA
<b>Reporting units</b>	ng/L	ng/L
<b>Reporting comments (interpretation)</b>	In an acute setting, Heart Failure is unlikely if NT-proBNP <300 ng/L Heart Failure is likely if: NT-proBNP >450 ng/L for patients <50 years of age NT-proBNP >900 ng/L for patients 50-75 years of age NT-proBNP >1800 ng/L for patients >75 years of age 2017 CCS HF Guidelines, CJC 2017	In an acute setting in the presence of appropriate clinical evaluation, the diagnosis of heart failure is: BNP < 100 ng/L: Unlikely BNP 100-400 ng/L: Possible, but other diagnoses must be considered. BNP > 400 ng/L: Very likely. 2017 CCS HF Guidelines, CJC 2017



**Appendix 1: 2-Hour Chest Pain Pathway for Quidel® TriageTrue® High Sensitivity Troponin-I Assay**





**Note:**

\*For all patients with abnormal hs-cTnI results, check the medical record for prior results. Many patients have stable abnormalities in hs-cTnI and measured concentrations similar to the patient’s baseline are reassuring.

Per European Society of Cardiology (ESC) Guidelines and 4<sup>th</sup> universal definition of MI, if the patient is >6h from symptom onset, has a hs-cTnI <21 ng/L (99<sup>th</sup> percentile upper reference limit), are pain-free, and have a low-risk presentation, they can be considered ruled out.

However, coronary ischemia has not been definitively excluded and unstable angina must be considered. Disposition after a single hs-cTnI <21 ng/L should only be considered for patients with low-risk clinical presentations who are >6 hours since symptoms onset, and should be used cautiously.

Troponin concentration may be elevated in the presence of kidney dysfunction. Patients with an eGFR <60 and an elevated troponin concentration should undergo serial testing to confirm a rising troponin concentration consistent with myocardial injury. Patients with an eGFR <60 who have troponin concentrations in the rule-out zone can be considered to have myocardial injury safely excluded.

<b>HEART Score Calculation</b>			
<b>History</b>	Highly suspicious		<b>2</b>
	Moderately suspicious		<b>1</b>
	Slightly suspicious		<b>0</b>
<b>ECG</b>	Significant ST-depression		<b>2</b>
	Non-specific repolarization disturbance, LBBB, LVH, Paced		<b>1</b>
	Normal		<b>0</b>
<b>Age</b>	≥ 65 years		<b>2</b>
	45-64 years		<b>1</b>
	≤ 44 years		<b>0</b>
<b>Risk Factors</b>	<input type="checkbox"/> Diabetes <input type="checkbox"/> Current smoker <input type="checkbox"/> Obesity <input type="checkbox"/> Family hx CAD <input type="checkbox"/> HTN (diagnosed) <input type="checkbox"/> HL (diagnosed)	≥ 3 risk factors or history of atherosclerotic disease	<b>2</b>
		1 or 2 risk factors	<b>1</b>
		No risk factors known	<b>0</b>
<b>hs-cTnI (peak)</b>	> 3x normal limit (64 ng/L or greater)		<b>2</b>
	1-3x normal limit (21-63 ng/L)		<b>1</b>
	< normal limit (< 21 ng/L)		<b>0</b>
<b>Total (10 maximum)</b>			
<b>HEART Score Interpretation</b>			
<b>Low Risk</b>			<b>0-3</b>
<b>Moderate Risk</b>			<b>4-6</b>
<b>High Risk</b>			<b>7-10</b>