











Leaders in Laboratory Medicine

- When only an Anatomic Pathology (AP) biopsy or Non-BAL Cytopathology specimen is collected (no BAL specimens), do not utilize the *Bronchoscopy Specimen Destination Label*. AP specimen will be tested at the collection site AP laboratory.
- When an Anatomic Pathology (AP) biopsy or Non-BAL Cytopathology specimen is collected in conjunction with BAL specimens follow the packaging instructions as outlined below. Use the *Bronchoscopy Specimen Destination Label*.
- When only a BAL Cytopathology specimen is collected, follow the packing instructions as outlined below. Use the *Bronchoscopy Specimen Destination Label*.

Order	Microbiology, Cytopathology, Hematology (ILD DIFF), Flow Cytometry, ProvLab	Anatomic Pathology Tissue Biopsy Specimens	Cytopathology Non-BAL specimens (e.g. TBNA)
Requisition	Bronchoalveolar Lavage/Bronchial Wash REQ9024MI  <ul style="list-style-type: none"> • Single or multiple specimens from one body site. Each different body site requires its own requisition. Do not use one requisition for multiple body sites. • All tests can be ordered on this requisition with the exception of: <ul style="list-style-type: none"> ➢ Tissue Biopsy specimens ➢ Non-BAL Cytopathology specimens • Include a contact name and number to resolve any specimen labelling or requisition problems. 	Anatomic Pathology Biopsy REQ19027 	Cytopathology Non-BAL specimen REQ9041CY-NON 
Sample Container	  		
Packaging	Ensure specimen container lids are secure to prevent any leakage. Place each specimen into an individual resealable specimen transport bag. Place requisition in outer pouch of individual resealable specimen transport bag.		
Prepare for transport to laboratory	Place all the individual specimen transport bags from same patient and specimen body site into a large resealable plastic bag labelled with the <i>Bronchoscopy Specimen Destination Label</i> REG5285 . DO NOT package specimens from multiple patients in the same large transport bag. Once packaged, immediately transport specimens to laboratory. DO NOT refrigerate. DO NOT transport bronchoscopy specimens using the Pneumatic Tube System.		