



**Malarial Smear Preparation**

**Sample Collection**      **5 mL EDTA Blood (optimal for adults).** Minimum sample volume for adults is 2 mL. Blood smears should be made as soon as possible after the blood is collected and ideally no later than 2 hours after collection. EDTA causes red blood cells to enlarge and lyse and will distort parasites making speciation difficult.

**Smears to Make**

- **6 thin smears**
- **4 thick smears**

**Materials**


- Semi-automatic pipetting device (e.g. Eppendorf) with ability to deliver 12 ul
- Clean frosted slides
- Micro hematocrit tube, disposable pipette, Diff SAFE/Smear Safe

**THIN Smears**      Make the same as regular Hematology smears. A well-prepared film is thick at one end and thin at the other (one layer of evenly distributed RBCs with no cell overlap).

Step	Action
1.	Mix EDTA blood well.
2.	Using a micro hematocrit tube, Diff Safe or disposable pipette, place a small drop of blood (approximately 2 – 3 mm diameter) about 1 cm from the frosted end of the slide.
3.	Using another slide spread the blood by drawing the slide back until it touches the drop of blood.
4.	Holding the spreader at approximately a 30-45° angle, smoothly, evenly draw it back into the drop of blood until the blood spreads about ¾ across the spreaders edge. Then gently pressing down on the spreader slide, push away from the drop with a smooth even motion, pushing the blood out along the length of the slide The finished smear should be at least 2.5 cm long with the tail ending approximately 1 cm from the end of the slide and have a monolayer of cells in 1/3 to 1/4 of end of slide with free margins on each side. If smears appear too thick, decrease the angle at which the spreader is held.
5.	Let dry 10 to 15 minutes prior packaging up.

**THICK Smears**

Step	Action
1.	Mix EDTA blood well
2.	Using a semi-automated pipette, dispense 12 ul blood onto the center of the slide. The use of 12 ul blood volume has been proven as optimal; when a larger volume drop is used there is a potential for the blood to wash off in the staining process. If a semi-automated pipette is not available use a Diff Safe, micro hematocrit tube or disposable pipette to deliver 2 small drops (about twice as much needed for a thin blood film).
3.	Spread blood out to the size of approximately a dime (15 to 20 mm in diameter) by using a stick or corner of another slide and spread the blood in a circular motion starting from the centre of the drop. The smear should be uniform in thickness and small newsprint should be just visible through the blood film.  If the slide does not appear to be of uniform thickness, gently rock slide to distribute blood evenly.
4.	Let dry 30 minutes prior to packaging up.



**SAMPLE DELIVERY**      Send all smears AND EDTA blood by STAT transportation to Calgary Microbiology at the Diagnostic and Scientific Centre (DSC).

**Questions? Contact Calgary Microbiology at the Diagnostic and Scientific Centre (DSC) 403-770-3278**