

This document is applicable at site(s):

**All Sites**

**Applicability** This document applies to all personnel of AHS Laboratory Services, the Lamont Health Centre and laboratories administered by Covenant Health.

**Purpose** This document provides instruction on special coagulation aliquoting.

**Procedure** **A. Aliquoting Chart**

**NOTE: A requisition must accompany the sample for all test codes indicated with an asterisks (\*)**

Test Code	Test Description	Acceptable Specimen Tube Reject all other specimen types	Frozen Aliquots Plasma = citrated
*APA	Antiphospholipid Antibodies or Cardiolipin Antibodies	<ul style="list-style-type: none"> <li>• Gold top tube</li> <li>• Ordered in conjunction with LUP.</li> </ul>	2 x 0.5 ml <b>serum</b> Mark <b>serum</b> on label
*APCE	APC Resistance or Activated Protein C Resistance or clot based Factor V Leiden	<ul style="list-style-type: none"> <li>• Glass/Plastic blue top tube</li> </ul>	2 x 1.0 ml plasma
APCGN and/or PROM	Factor V Leiden mutation and/or Prothrombin 20210 (or Factor II) mutation	<ul style="list-style-type: none"> <li>• 1 – 4 mL mauve top tube @ 4°C</li> <li>• MUST be sterile, dedicated sample</li> </ul>	N/A – MUST be submitted as whole blood. Send Cold.
COAGZ	Special Coagulation Evaluation: This is a Morphology/Coagulation Pathologist On-Call initiated order - they will contact UAH Special Coag with specific test order	<ul style="list-style-type: none"> <li>• Glass/Plastic blue top tube</li> </ul>	3 x 0.5 ml <b>Platelet Free Plasma</b> min. Mark <b>PFP</b> on label
*F2	Factor II	<ul style="list-style-type: none"> <li>• Glass/Plastic blue top tube</li> </ul>	2 x 0.5 ml plasma per test
*F5	Factor V		
*F7	Factor VII		
*F8	Factor VIII		
*F8INH	Factor VIII Inhibitor	<ul style="list-style-type: none"> <li>• Glass/Plastic blue top tube</li> </ul>	2 x 1.0 ml plasma
*RISTO	Ristocetin Cofactor	<ul style="list-style-type: none"> <li>• Glass/Plastic blue top tube</li> </ul>	2 x 0.5 ml plasma per test
*VWFAG	Von Willebrand Factor <b>or</b> VWF Antigen		
*F9	Factor IX	<ul style="list-style-type: none"> <li>• Glass/Plastic blue top tube</li> </ul>	2 x 0.5 ml plasma per test
*F10	Factor X		
*F11	Factor XI		
*F12	Factor XII		
*F13	Factor XIII	<ul style="list-style-type: none"> <li>• Glass/Plastic blue top tube</li> </ul>	2 x 1.0 ml plasma
*HITB	Heparin induced thrombocytopenia assay	<ul style="list-style-type: none"> <li>• Gold top tube</li> </ul>	2 x 0.5 ml <b>serum</b> Mark <b>serum</b> on label
*LUP	Lupus anticoagulant <b>or</b> Dilute Russell's Viper Venom Time	<ul style="list-style-type: none"> <li>• 2 GLASS blue top tubes</li> <li>• <b>Plastic WILL NOT be accepted</b></li> <li>• Ordered in conjunction with APA</li> </ul>	2 x 1.0 ml <b>Platelet Free Plasma</b> Mark <b>PFP</b> on label
HEPRN	Low Molecular Weight Heparin <b>or</b> Anti-Xa assay <b>or</b> Unfractionated Heparin (LMWH: Enoxaparin, Fragmin, Dalteparin)	<ul style="list-style-type: none"> <li>• Glass/Plastic blue top tube</li> <li>– <b>must be collected 4-6 hours post dose. Room Temp.</b></li> </ul>	MUST be processed within 4 hrs of collection. 2 x 0.5 ml <b>Platelet Free Plasma</b> Mark <b>PFP</b> on label
AT3	Anti-thrombin	<ul style="list-style-type: none"> <li>• 1 Glass/Plastic blue top tube</li> </ul>	1 x 0.5 mL plasma
*PROTC	Protein C	<ul style="list-style-type: none"> <li>• 2 Glass/Plastic blue top tubes</li> </ul>	2 x 0.5 ml plasma per test
*PROTS	Protein S		
*PTIN	PT Inhibitor Screen	<ul style="list-style-type: none"> <li>• 2 Glass/Plastic blue top tubes</li> </ul>	2 x 2.0 ml plasma
*PTTIN	PTT Inhibitor Screen	<ul style="list-style-type: none"> <li>• 3 Glass/Plastic blue top tubes</li> </ul>	2 x 3.0 ml <b>Platelet Free Plasma</b> Mark <b>PFP</b> on label
TT	Thrombin Time	<ul style="list-style-type: none"> <li>• 1 Glass/Plastic blue top tube</li> </ul>	2 x 0.5 ml plasma

**B. Aliquoting Instructions**

Step	Detail	Information																				
1. Spin specimen.	1.1) 2500 g/rcf for 10 minutes. 1.2) If specimen received spun proceed to <b>Step 2</b> . 1.3) <b>Aliquot and freeze sample aliquots within 4 hrs of collection.</b>	<ul style="list-style-type: none"> <li>Platelet poor plasma must have a platelet count of <math>&lt;10 \times 10^9/L</math>.</li> </ul>																				
2. Check specimen integrity before pooling.	2.1) Ensure patient identification for all collection tubes are the same. 2.2) Using a <b>plastic</b> pipette, remove the plasma or serum from cells one tube at a time – avoid the platelet and WBC layer when removing plasma and/or clot separator when removing serum. <table border="1" data-bbox="284 577 1144 1050"> <thead> <tr> <th>If...</th> <th>Then...</th> </tr> </thead> <tbody> <tr> <td>Removing plasma</td> <td> <ul style="list-style-type: none"> <li>Keep plasma in pipette.</li> <li>Check sample for a clot prior to adding to the pool.</li> </ul> <table border="1" data-bbox="446 672 1128 913"> <thead> <tr> <th>If...</th> <th>Then...</th> </tr> </thead> <tbody> <tr> <td>Clotted</td> <td> <ul style="list-style-type: none"> <li>DO NOT add to the plasma pool.</li> <li>Discard plasma.</li> </ul> </td> </tr> <tr> <td>Hemolyzed</td> <td> <ul style="list-style-type: none"> <li>Freeze hemolyzed plasma separately and label as HEMO.</li> </ul> </td> </tr> <tr> <td>Not clotted and not hemolyzed</td> <td> <ul style="list-style-type: none"> <li>Proceed to <b>Step 3</b>.</li> </ul> </td> </tr> </tbody> </table> </td> </tr> <tr> <td>Removing serum</td> <td> <table border="1" data-bbox="446 913 1128 1050"> <thead> <tr> <th>If...</th> <th>Then...</th> </tr> </thead> <tbody> <tr> <td>Hemolyzed</td> <td> <ul style="list-style-type: none"> <li>Freeze hemolyzed plasma separately and label as HEMO.</li> </ul> </td> </tr> <tr> <td>Not hemolyzed</td> <td> <ul style="list-style-type: none"> <li>Proceed to <b>Step 3</b>.</li> </ul> </td> </tr> </tbody> </table> </td> </tr> </tbody> </table>	If...	Then...	Removing plasma	<ul style="list-style-type: none"> <li>Keep plasma in pipette.</li> <li>Check sample for a clot prior to adding to the pool.</li> </ul> <table border="1" data-bbox="446 672 1128 913"> <thead> <tr> <th>If...</th> <th>Then...</th> </tr> </thead> <tbody> <tr> <td>Clotted</td> <td> <ul style="list-style-type: none"> <li>DO NOT add to the plasma pool.</li> <li>Discard plasma.</li> </ul> </td> </tr> <tr> <td>Hemolyzed</td> <td> <ul style="list-style-type: none"> <li>Freeze hemolyzed plasma separately and label as HEMO.</li> </ul> </td> </tr> <tr> <td>Not clotted and not hemolyzed</td> <td> <ul style="list-style-type: none"> <li>Proceed to <b>Step 3</b>.</li> </ul> </td> </tr> </tbody> </table>	If...	Then...	Clotted	<ul style="list-style-type: none"> <li>DO NOT add to the plasma pool.</li> <li>Discard plasma.</li> </ul>	Hemolyzed	<ul style="list-style-type: none"> <li>Freeze hemolyzed plasma separately and label as HEMO.</li> </ul>	Not clotted and not hemolyzed	<ul style="list-style-type: none"> <li>Proceed to <b>Step 3</b>.</li> </ul>	Removing serum	<table border="1" data-bbox="446 913 1128 1050"> <thead> <tr> <th>If...</th> <th>Then...</th> </tr> </thead> <tbody> <tr> <td>Hemolyzed</td> <td> <ul style="list-style-type: none"> <li>Freeze hemolyzed plasma separately and label as HEMO.</li> </ul> </td> </tr> <tr> <td>Not hemolyzed</td> <td> <ul style="list-style-type: none"> <li>Proceed to <b>Step 3</b>.</li> </ul> </td> </tr> </tbody> </table>	If...	Then...	Hemolyzed	<ul style="list-style-type: none"> <li>Freeze hemolyzed plasma separately and label as HEMO.</li> </ul>	Not hemolyzed	<ul style="list-style-type: none"> <li>Proceed to <b>Step 3</b>.</li> </ul>	<b>NOTE:</b> Check for a clot by rimming the red cell button with an applicator stick.
If...	Then...																					
Removing plasma	<ul style="list-style-type: none"> <li>Keep plasma in pipette.</li> <li>Check sample for a clot prior to adding to the pool.</li> </ul> <table border="1" data-bbox="446 672 1128 913"> <thead> <tr> <th>If...</th> <th>Then...</th> </tr> </thead> <tbody> <tr> <td>Clotted</td> <td> <ul style="list-style-type: none"> <li>DO NOT add to the plasma pool.</li> <li>Discard plasma.</li> </ul> </td> </tr> <tr> <td>Hemolyzed</td> <td> <ul style="list-style-type: none"> <li>Freeze hemolyzed plasma separately and label as HEMO.</li> </ul> </td> </tr> <tr> <td>Not clotted and not hemolyzed</td> <td> <ul style="list-style-type: none"> <li>Proceed to <b>Step 3</b>.</li> </ul> </td> </tr> </tbody> </table>	If...	Then...	Clotted	<ul style="list-style-type: none"> <li>DO NOT add to the plasma pool.</li> <li>Discard plasma.</li> </ul>	Hemolyzed	<ul style="list-style-type: none"> <li>Freeze hemolyzed plasma separately and label as HEMO.</li> </ul>	Not clotted and not hemolyzed	<ul style="list-style-type: none"> <li>Proceed to <b>Step 3</b>.</li> </ul>													
If...	Then...																					
Clotted	<ul style="list-style-type: none"> <li>DO NOT add to the plasma pool.</li> <li>Discard plasma.</li> </ul>																					
Hemolyzed	<ul style="list-style-type: none"> <li>Freeze hemolyzed plasma separately and label as HEMO.</li> </ul>																					
Not clotted and not hemolyzed	<ul style="list-style-type: none"> <li>Proceed to <b>Step 3</b>.</li> </ul>																					
Removing serum	<table border="1" data-bbox="446 913 1128 1050"> <thead> <tr> <th>If...</th> <th>Then...</th> </tr> </thead> <tbody> <tr> <td>Hemolyzed</td> <td> <ul style="list-style-type: none"> <li>Freeze hemolyzed plasma separately and label as HEMO.</li> </ul> </td> </tr> <tr> <td>Not hemolyzed</td> <td> <ul style="list-style-type: none"> <li>Proceed to <b>Step 3</b>.</li> </ul> </td> </tr> </tbody> </table>	If...	Then...	Hemolyzed	<ul style="list-style-type: none"> <li>Freeze hemolyzed plasma separately and label as HEMO.</li> </ul>	Not hemolyzed	<ul style="list-style-type: none"> <li>Proceed to <b>Step 3</b>.</li> </ul>															
If...	Then...																					
Hemolyzed	<ul style="list-style-type: none"> <li>Freeze hemolyzed plasma separately and label as HEMO.</li> </ul>																					
Not hemolyzed	<ul style="list-style-type: none"> <li>Proceed to <b>Step 3</b>.</li> </ul>																					
3. Pool plasma and/or serum.	3.1) Into plastic container.	<ul style="list-style-type: none"> <li>Accuvette or sterile urine container.</li> </ul>																				
4. Transfer the pooled plasma.	4.1) Mix the pool, aliquot into labeled plastic tubes, and freeze at -70°C – refer to <b>Section A</b> . <b>** Preferred tube:</b> 4 mL VWR #718585 and screw top lid VWR #718593 or refer to <i>Client Supply Order – Edmonton Zone Form**</i>	<b>NOTE:</b> Sample may be stored at -20°C for a short period (eg. 2 days) in a NON-frost free freezer; sample should be shipped on dry ice.																				

**C. Preparation of Platelet Free Plasma (PFP)**

Step	Detail	Information
1. Perform the primary spin.	1.1) Follow steps 1-3 in <b>Section B</b> .	
2. Perform a secondary spin.	2.1) Mix the pool and aliquot into labeled plastic eppendorf tubes. 2.2) Spin in a high speed centrifuge for 6 mins. <b>NOTE: If a high speed centrifuge is not available, transfer plasma to plastic tubes and repeat Step 1 in Section B.</b>	<ul style="list-style-type: none"> <li>UAH (Eppendorf 5415D – 12000 rpm /13600 g/rcf – 6 min)</li> <li>RAH (Eppendorf 5452 – 13400 rpm – 6 min)</li> </ul>
3. Label aliquot tubes.	3.1) Indicate on the label that it is <b>PFP</b> on appropriate number of tubes. <b>NOTE: Aliquot tubes must be plastic – refer to Section B, Step 4.</b>	
4. Aliquot PFP.	4.1) Remove plasma with <b>plastic</b> pipette and pool – be careful not to disturb the cell button. 4.2) Mix the pool, aliquot into labeled plastic tubes, and freeze at -70°C – refer to <b>Section A</b> .	<b>NOTE:</b> Sample may be stored at -20°C for a short period (ie. 2 days) in a NON-frost free freezer; sample should be shipped on dry ice.

**References** N/A

**Related Documents**

Current Version of:

Document	Document Control Number
Client Supply Order – Edmonton Zone Form	RRRSOR00005

**Contact Information**

University of Alberta Hospital Routine Coagulation (780) 407-7232