

Leaders in Laboratory Medicine

Laboratory Bulletin

| DATE: | 2022 April 4 |
|-------|--|
| TO: | Peace Region Hospitals that Refer to Grande Prairie Regional Hospital |
| FROM: | A. Brent Mendez, Regional Lab Medicine Site Chief, GPRH, Grande Prairie, Alberta |
| RE: | Body Fluid Tube Type Change |

PLEASE POST OR DISTRIBUTE AS WIDELY AS APPROPRIATE

Key Message

 All body fluids submitted for chemical analysis to the Grande Prairie Regional Hospital (GPRH) need to be placed in a sterile container, although red top (no additive) tubes and gold top (serum separator tube) tubes are also acceptable.

Background

• Current acceptable tube type is dark green (sodium heparin). In order to provinicially align practices and following biochemistry recommendations, our acceptability criteria is changing to a sterile container, however, a red top or gold top tube is also acceptable.

How this will impact you

Be aware of the changing acceptability requirements.

Effective immediately

Questions/Concerns

- Dustin Proctor, PhD, Clinical Biochemist, Alberta Precision Laboratories.
- Allen Brent Mendez, M.D., Alberta Precision Laboratories, Regional Lab Medicine Site Chief, GPRH, Grande Prairie, Alberta

Approved by

• Kareena Schnabl, MSc, PhD, FCACB, Section Chief Biochemistry North Sector, Clinical Biochemist, Alberta Precision Laboratories.