



Laboratory Bulletin

Lead	lers	in l	Lab	ora	tory	Me	dic	ine
Loud	1010			, Oi a	LOI y	IVIC	aic	1110

DATE:	2022 March 21
TO:	All Healthcare Professionals in Central, Edmonton, and North Zone
FROM:	Alberta Precision Laboratories (APL) North Sector, and DynaLIFE Medical Labs
RE:	Changes and Updates to Endocrine Tests

PLEASE POST OR DISTRIBUTE AS WIDELY AS APPROPRIATE

Key Message

• Effective **March 29, 2022**, there will be various ordering or reporting changes for the following tests: Aldosterone; Aldosterone/renin ratio (ARR); Aldosterone, urine, 24 hour; Thyroglobulin panel (thyroglobulin and anti-thyroglobulin); and Anti-thyroid peroxidase (TPO).

Background

- These changes will help optimize test interpretation and improve ordering practices.
- They are also important for supporting laboratory standardization efforts across the province.

How this will impact you

Test	Change details	Performing site
Aldosterone	 Preferred blood collection container type will switch from serum separator tube (gold top) to EDTA plasma (lavender top). This will improve reliability of high aldosterone concentrations, consolidate collection of aldosterone and renin into a single container, and prevent delays in ARR reporting. 	University of Alberta Hospital (UAH)
ARR	 Aldosterone >400 pmol/L will be removed as criteria for calculating and reporting the ARR. This change will ensure that a subset of patients with aldosterone ≤400 pmol/L will not be missed for primary hyperaldosteronism screening. 	UAH
Aldosterone, urine, 24 hour	 The reference interval for adults will append as a comment. A comment will display for patients <18 years indicating that no reference interval is available. 	UAH
Thyroglobulin panel	 Thyroglobulin and anti-thyroglobulin will be combined into a single panel. Neither test will be available as an individual order. This change will allow all thyroglobulin results to be assessed for the potential of anti-thyroglobulin interference in the same report. 	UAH
TPO	 Measuring units will change from kU/L to the SI units kIU/L. Numerical values are equivalent and no conversion is required. The reference interval will change from <50 kIU/L to <35 kIU/L to better reflect current methodology. The upper reportable limit will be 600 kIU/L. Results higher than 600 kIU/L will be reported as >600 kIU/L. 	UAH and DynaLIFE Medical Labs



Leaders in Laboratory Medicine

Action Required

- Be aware of test updates and reporting changes as outlined in the above table.
- Laboratory Test Directory & Collection Information can be accessed at https://www.albertahealthservices.ca/lab/Page3217.aspx

Questions/Concerns

- Dr. Josh Raizman, Clinical Biochemist, APL, (780) 718-2402, josh.raizman@aplabs.ca
- Dr. Albert Tsui, Clinical Biochemist, APL, (587) 782-2674, albert.tsui@aplabs.ca
- DynaLIFE Clinical Chemists, (780) 451-3702 ext. 3572, clinical.chemists@dynalife.ca

Approved by

- Dr. Kareena Schnabl, Section Chief, Clinical Biochemistry, APL North Sector
- Dr. Michael Mengel, Medical Director, APL North Sector
- Dr. Mathew Estey, Director of Clinical Chemistry, DynaLIFE Medical Labs
- Dr. Erene Farag, Interim Medical Director, DynaLIFE Medical Labs