

FOR THE DONNING PROCEDURE USE A CHALLENGE RESPOND METHOD

CHALLENGE RESPOND = HCW completes checklist actions.
 After completion of sequence, **Buddy** reads off checklist item (CHALLENGE)
HCW confirms they have completed checklist item (RESPOND)

1. PREPARATION <ul style="list-style-type: none"> • Pull hair back and secure it off neck and face • Remove all extra items on person 	<input type="checkbox"/> <input type="checkbox"/>
2. HAND HYGIENE <ul style="list-style-type: none"> • Perform hand hygiene using alcohol based hand rub If your hands look or feel dirty use soap and water.	<input type="checkbox"/>
3. COVERALLS TO SHOULDERS <ul style="list-style-type: none"> • Step into coveralls, pull up to shoulders • Put arms in arm holes Do not put on hood, do not zip up.	<input type="checkbox"/> <input type="checkbox"/>
4. SHOE COVERS <ul style="list-style-type: none"> • Pull shoe covers over feet of coveralls 	<input type="checkbox"/>
5. HAND HYGIENE <ul style="list-style-type: none"> • Perform hand hygiene, using alcohol based hand rub 	<input type="checkbox"/>
6. INNER GLOVES (SURGICAL) <ul style="list-style-type: none"> • Put on surgical gloves. Cuffs of gloves go under cuffs of coveralls 	<input type="checkbox"/>
7. MASK OR N95 RESPIRATOR <ul style="list-style-type: none"> • Remove personal eye glasses, if worn 7A. MASK <ul style="list-style-type: none"> • Put on mask • Put ear loops over ears • Pinch nose piece to fit bridge of nose using both hands 7B. N95 RESPIRATOR <ul style="list-style-type: none"> • Stretch out elastic bands; cup N95 in one hand and hold to face • Put top elastic on crown of head put bottom elastic on base of head • Ensure bands are not crossed; no hair is under edge of mask • Conduct a seal check <ul style="list-style-type: none"> • Put on personal eye glasses, or optional goggles 	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
8. BOUFFANT HEAD COVER <ul style="list-style-type: none"> • Put on bouffant head cover • Ensure ears and all hair is covered (use mirror or another person) 	<input type="checkbox"/> <input type="checkbox"/>

Coverall Donning Checklist
 Viral Hemorrhagic Fever (Ebola)
 Infection Prevention and Control - September 2015

<p>9. COVERALLS-HOOD AND ZIPPER SEAL</p> <ul style="list-style-type: none"> • Pull up coverall hood • Fully zip up coveralls • Buddy peels backing off coverall flap • Buddy holds flap out and straight • HCW presses flap closed, ensuring there are no folds or bubbles • HCW or Buddy press firmly on the top of the flap for a few seconds to seal it well 	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<p>10. FULL-NECK HOOD</p> <ul style="list-style-type: none"> • Pull hood over head • Buddy ensures neck and bouffant are covered and visual field is clear 	<input type="checkbox"/> <input type="checkbox"/>
<p>11. APRON</p> <ul style="list-style-type: none"> • Put on apron. Make sure apron covers from neck to knees and to wrists • Insert thumbs in loops of sleeves • Tie apron at back of waist • Buddy ensures apron covers full neck hood at shoulders 	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<p>12. FACE SHIELD</p> <ul style="list-style-type: none"> • Put on face shield • Adjust to fit square on face 	<input type="checkbox"/> <input type="checkbox"/>
<p>13. OUTER GLOVES</p> <ul style="list-style-type: none"> • Put on nitrile or a second pair of surgical gloves over cuffs of apron • Buddy to assist if necessary 	<input type="checkbox"/> <input type="checkbox"/>
<p>14. INSPECTION</p> <ul style="list-style-type: none"> • Buddy inspects to ensure HCW PPE is secure and all closures are flat. Before HCW enters room, write HCW name and occupation, and time on front and arm of apron. <p>To ensure PPE stays in place before entering room, HCW should:</p> <ul style="list-style-type: none"> • bend forward and gently shake head • perform a squat • fully extend arms overhead 	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<p>PERFORM A DOUBLE CHECK TO ENSURE DONNING OF PPE IS COMPLETE</p>	