

Using A Voice Amplifier

Description

A voice amplifier is a portable device that helps others to hear a person's speech. The person speaks into a microphone that amplified the sound so that others can hear what they have said. Most units can be used with any type of external microphone. There are a variety of device options to consider depending on the person's needs, including Bluetooth VoiceSaver[®] ; ChatterVox[®] ProH; Voice Buddy[®] Mini Buddy ; Spokesman[®] Voice Amplifier and JustAMP[®] Digital Amplifier.



Who May Benefit?

People with a quiet, weak voice, including those with respiratory insufficiency and those who fatigue easily from speaking.

How to Use

- Attach the speaker/ battery pack to the person's waist or wheelchair
- Turn on the device. Position the microphone and adjust the volume to fit the situation

How to Use

Ask what, if any support the person needs to care for and use the device:

- Turning it on and off
- Adjusting the volume
- Putting on and positioning the microphone
- Recharging or replacing the batteries on a regular basis
- Ensuring that the device and charger/spare batteries are stored in its case in a safe location to avoid loss or damage
- Making sure it is accessible and that family, friends and healthcare providers are encouraged to use it during interactions

