

Collection Instructions for: *Chlamydia Trachomatis* and *Neisseria Gonorrhoeae* NAAT Testing

For specimens from the **cervix, urethra or eye** use:

UNISEX Collection kit (BLUE swab; WHITE label)



For specimens from the **throat, vagina or rectum** use:

MULTITEST Collection kit (PINK swab; ORANGE label)



For **urine** use:

URINE collection kit (yellow label)



For testing of specimens from other sites, call the Microbiologist on Call to discuss collection.

Calgary – 403-944-1200

Edmonton – 780-407-7121