

Client Supply Order – Provincial Form

Fax Completed Form to: UAH Reagent Room Fax №: 780-407-8599 Phone №: 780-407-8401

Site/Address:	Contact Name:	Phone №:
	Order Date:	Fax №:

For testing/collection information contact Client Response: **Phone №:780-407-7484 Toll Free: 1-866-407-7767**

This form can be downloaded from: <http://www.albertahealthservices.ca/assets/wf/lab/wf-lab-client-supply-order-provincial-form.pdf>

NOTE: Reagent Room is open daily Monday to Friday. Two to three working day's notice is required to prepare and ship orders.

Description	Unit of Issue	Quantity	Shipped
Bone Marrow Media (BMM)	Each		
Catecholamine collection tube	Each		
Glucagon collection tubes	Each		
Green 4.5 mL Lithium Heparin collection tube	Each		
Green 4 mL Sodium Heparin collection tube	Each		
Grey 4 mL Sodium Fluoride Potassium Oxalate collection tube	Each		
Grey 6 mL Sodium Fluoride Potassium Oxalate collection tube	Each		
PAXgene Blood RNA collection tubes	Each		
PAXgene Bone Marrow collection tubes	Each		
Royal Blue EDTA 7 mL collection tubes	Each		
Royal Blue Plain 7 mL collection tubes includes Sarstedt Pour-off tube	Each		
Salivette collection tubes for Salivary Cortisol	Each		
Sarstedt Pour-off tubes for Trace Elements, Urine	Each		
Special Coagulation pour-off tubes with lids	Pkg of 100		
Stool Can, 72 hr., pre-weighed	Each		
TCA (Trichloroacetic Acid), 5%	50mL		
Tissue Transport Media (TTM)	Each		
Urine Container, 24 hr., Amber, with Acid	Each		
Yellow ACD/A 8.5 mL collection tube	Each		

UAH Reagent Room use Only

Date Received:	Date Shipped:	Number of Pieces:
-----------------------	----------------------	--------------------------