

IGFBP-3 Tanner Stage Reference Ranges

Effective Date: February 28, 2011

Chemistry Client Resource

Version: 1.0

 This document is applicable at site(s):

UAH

IGFBP-3 Tanner Stage Reference Ranges

IGFBP-3, mg/L

Tanner Stage	Median	Central 95% Range
--------------	--------	-------------------

 Female, $n = 226$

1	3.6	1.2-6.4
2	4.5	2.8-6.9
3	5.3	3.9-9.4
4	5.9	3.3-8.1
5	5.6	2.7-9.1

 Male, $n = 194$

1	3.6	1.4-5.2
2	3.9	2.3-6.3
3	5.4	3.1-8.9
4	6.5	3.7-8.7
5	5.2	2.6-8.6

Reference IGFBP-3 Kit for the Siemens Immulite 2000 package insert (PIL2KGB-12, 2008-07-29)